

# Synovetin OA

[Polysynovial Tin (<sup>117m</sup>Sn Colloid)]

Veterinary Device for Use in Dogs

## NAME Synovetin 30<sup>®</sup>

Sn (<sup>117m</sup>Sn) stannic colloid in suspension, pH 7.8, supplied as a 2-4 mL (2% - 140 MBq/g), suspension for intrasynovial (SI) injection.

## NET QUANTITY

This container is pre-filled with up to 4 mL (200 MBq) of the eluate and time to first use day.

1 mL of suspension contains 2-4 mL (24-140 MBq) of Sn (<sup>117m</sup>Sn) stannic colloid in suspension pH 7.8 for SI use and time of first use.

## PRODUCT DESCRIPTION

Synovetin 30<sup>®</sup> is a suspension sterile, isotonic veterinary device comprising a colloid stannic suspension with sulphuric acid and 0.9% w/v of NaCl. 95% of the particles have a size between 1.0 µm and 20 µm (D50) and water insoluble. The <sup>117m</sup>Sn emits monoenergetic conversion electrons (average energy: 137-150 keV) with a probability 100% and isotope gamma emission (511 keV) with 36% abundance, accompanying low energy emissions are Auger electrons (<10 keV) and X-rays (<20 keV). The half-life of <sup>117m</sup>Sn is 14 days. <sup>117m</sup>Sn decaying by isomeric transition emits <sup>117g</sup>Sn.

Synovetin 30<sup>®</sup> includes a suspension with stannic (<sup>117m</sup>Sn) stannic colloid (SI), conversion electrons (SI), and gamma radiation, and has no chemical.

## INDICATIONS FOR USE

Synovetin 30<sup>®</sup> is a veterinary device consisting of a homogeneous tin colloid with water soluble in 0.9% w/v low energy conversion electrons confined to the joint space. The colloid is composed of nanoparticles (1.0 µm to 20 µm) that are retained in the joint space of the dog. The particles are absorbed and retained by synovial and macrophages in the synovium, resulting in synthesis and retention of extracellular matrix. Stimulation of the pro-inflammatory cells release inflammatory cells in the joint space, thereby reducing pain associated with arthritis. The data, including radiographic evidence, supports use in Grade 1, 2, and 3 osteoarthritis (OA) of the elbow joint.

## CONTRA

Contra use includes this device to use by or on the order of a licensed veterinarian trained in the use of radioactive materials and products.

Use of this product is restricted to facilities with a compatible Radioactive Material (RM) license.

## INDICATIONS

Synovetin 30<sup>®</sup> is intended to reduce synovial and associated pain of canine elbow joint affected with osteoarthritis.

## WARNING

Do not exceed 4.0 mL (200 MBq) of relative activity per dog per treatment, but for use in humans. Keep this and all preparations out of reach of children. Consult authorities in case of accidental ingestion or injection by humans.

## PRECAUTIONS

Injection should be performed only by a licensed veterinarian skilled in the delivery of this article (SI) operators who is trained at a facility that has a RM license.

Operator should technique used for correct dosing injection.

## CONTRAINDICATIONS

Use the SI (SI) to determine the appropriate dose. These were determined using the elbow joint.

For example, a dog weighing 20 lbs. receives an SI dose of 0.2 mL (20 MBq) in each elbow to be treated.

Dog Weight (lbs.)	Synovetin 30 <sup>®</sup> Dose per Elbow Joint (mL)
10-15 lbs.	0.6
20-25 lbs.	0.8
30-35 lbs.	1.2
40-45 lbs.	1.5
50-55 lbs.	1.7
60-65 lbs.	1.8
70-75 lbs.	2.2
80-85 lbs.	2.4
90-95 lbs.	2.6
100-110 lbs.	2.8
110 lbs. and over	3.0

**These will be linked to 1.0 mL minimum joint elbow weight amounts 100 lbs., with the total body dose not exceeding 6.0 mL (600 MBq), two elbow joints in 10-15 lb. or greater sized dogs.**

## PREPARATION FOR USE

Synovetin 30<sup>®</sup> is provided in a 2 mL, glass vial with a lead cylinder. Each vial is for use with a syringe only.

The product should be stored in the cardboard shipping container until needed for use. The **prescribed dose** should be **administered on the date** noted on the container accompanying the Synovetin 30<sup>®</sup> because it can be administered the day before or after 1 consecutive month (up to 30 days) after use. Proper preparation for equipment and procedure for handling radioactive work products, including shielding, safety areas, safety program, leak-testing, gloves, lead aprons, booties, and respirator masks.

STEP 1. Allow vial to equilibrate the dose after 4 springs and prior to entering the sterile area around the lead cylinder (only **while the lead cylinder is approximately 90 seconds to ensure proper mixing** of the product).

STEP 2. Remove the vial cap from the lead cylinder and dispose of it appropriately.

STEP 3. Remove the lead cylinder (LC), but do not remove the glass vial from the lead cylinder.

STEP 4. Remove the vial from the cap from the LC and retain for placement on the dog after the dose is finished.

STEP 5. Attach a sterile syringe (2 mL, or other appropriate volume) to a 22-ga. needle. When placed, use a syringe (used to maintain sterility) withdraw dose so the vial is approximately 2/3 full and the vial is held in a sterile position.

STEP 6. While holding the container at an approximately 45° angle, insert the needle through the septum.